

Atomic Spectra

Lab: 07

Jaden Moore

Orange Coast College
Physics A285L
December 6th, 2021

1 Introduction

In this lab, we analyze the light emitted by sodium and hydrogen as it enters through a diffraction grating. We can then use various tools and diffraction grating techniques to obtain an experimental wavelength for the light and compare it to the theoretical values to obtain a percent error between the two.

2 Measuring the wavelength of Sodium Light Source

Given a diffraction grating of 600 slits per mm we can calculate the distance between each slit such that

$$d = \frac{1}{\text{slits per mm}} = \frac{1}{600} = 1667 \text{ nm} \quad (1)$$

From the experiment, we see two yellow lines a distance away from the central bright maximum. We can then measure the two angles between the yellow lines and the central bright maximum. From this, we obtain the wavelength of the incoming light using the following formula,

$$\lambda = \frac{d \sin \theta}{m} \quad (2)$$

Below we put into a table the measured angles from the 1st order bright spot where $m = 1$ and then calculate the experimental wavelength for the light using Equation (2).

Table 1: Sodium wavelength results for $m = 1$	
Angle θ [deg]	Wavelength λ [nm]
20.71	589.514
20.69	588.970

Using Table (1), we can then calculate the doublet width $\Delta\lambda$ such that

$$\Delta\lambda = |\lambda_1 - \lambda_2| = 589.514 - 588.970 = 0.544 \text{ nm} \quad (3)$$

From the lab, we get that the difference in energy between the two doublets is calculated from the formula

$$\Delta E = h(v_2 - v_1) = hc \left(\frac{1}{\lambda_2} - \frac{1}{\lambda_1} \right) \quad (4)$$

Where h is Planck's constant and c is the speed of the light. Using Equation (3) we get the difference in energy to be

$$\Delta E = (4.136 \times 10^{-15} \text{ eVs})(3.0 \times 10^8 \text{ m/s}) \left(\frac{1}{5.8897 \times 10^{-7} \text{ m}} - \frac{1}{5.8951 \times 10^{-7} \text{ m}} \right) = 1.930 \times 10^{-3} \text{ eV} \quad (5)$$

We can then calculate the magnetic field seen by the electron. That is,

$$B = \frac{\Delta E}{g\mu_b} = \frac{1.930 \times 10^{-3}}{2(5.798 \times 10^{-5})} = 16.64 \text{ T} \quad (6)$$

In this case, we notice that the magnetic field generated is 16.64T which is much stronger than that of a magnetic resonance imaging (MRI) machine which is around 1 T.

3 Measuring the Spectrum of Hydrogen

Now we consider the spectra produced by a hydrogen atom. Using the angles measured from the colored bands to the central bright maximum, we can calculate the wavelength of the light using Equation (2) and then obtain a percent error between the experimental and theoretical values. We can calculate the percent error between the two wavelengths such that

$$\% \text{ error} = \left(\frac{\lambda_{exp} - \lambda_{theory}}{\lambda_{theory}} \right) 100 \quad (7)$$

Below we put into a table the results collected from the experiment conducted on the hydrogen atom.

Table 1: Hydrogen wavelength results for $m = 1$				
Color	Angle θ [deg]	λ_{exp} [nm]	λ_{theory} [nm]	% error
Red	23.10	654.025	656.3	0.346
Aqua	16.95	485.992	486.1	0.022
Blue	15.09	433.980	434.0	0.004
Violet	14.24	410.056	410.2	0.035
Ultra Violet	13.77	396.787	397.0	0.054

In this case, we see that the percent error is within 1% of the theoretical values.

4 Conclusion

From this lab, we get that the light emitted from sodium and hydrogen produces spectral lines when incident on a diffraction grating. Thus, we can calculate the wavelength of the light emitted by the elements by considering the formulas associated with diffraction grating. We found that the wavelengths obtained from the experiment were within 1% error of the theoretical values which indicate that the experiment was a success. The biggest takeaway of the lab is gaining a better understanding of how we can use methods like diffraction grating to measure light, even at an atomic scale.